



SV2B Antibody, FITC conjugated

Product Code	CSB-PA022979LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L1I2
Immunogen	Recombinant Human Synaptic vesicle glycoprotein 2B protein (1-108AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptic vesicle glycoprotein 2B, SV2B, KIAA0735
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SV2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.