



SWAP70 Antibody, Biotin conjugated

Product Code	CSB-PA887026LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UH65
Immunogen	Recombinant Human Switch-associated protein 70 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Switch-associated protein 70 (SWAP-70), SWAP70, KIAA0640
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SWAP70
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.