



TAGAP Antibody, Biotin conjugated

Product Code	CSB-PA818676LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N103
Immunogen	Recombinant Human T-cell activation Rho GTPase-activating protein (97-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell activation Rho GTPase-activating protein (T-cell activation GTPase-activating protein), TAGAP, TAGAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TAGAP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.