



TARS2 Antibody, Biotin conjugated

Product Code	CSB-PA023132LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW92
Immunogen	Recombinant Human Threonine--tRNA ligase, mitochondrial protein (369-718AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Threonine--tRNA ligase, mitochondrial (EC 6.1.1.3) (Threonyl-tRNA synthetase) (ThrRS) (Threonyl-tRNA synthetase-like 1), TARS2, TARSL1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TARS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.