



TBC1D20 Antibody, Biotin conjugated

Product Code	CSB-PA846594LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BZ9
Immunogen	Recombinant Human TBC1 domain family member 20 protein (1-237AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 20, TBC1D20, C20orf140
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	TBC1D20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.