



# TBC1D9B Antibody, HRP conjugated

Product Code	CSB-PA724138LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q66K14
Immunogen	Recombinant Human TBC1 domain family member 9B protein (724-1023AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 9B, TBC1D9B, KIAA0676
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D9B