




# TCEAL7 Antibody

<b>Product Code</b>	CSB-PA863109LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BRU2
<b>Immunogen</b>	Recombinant Human Transcription elongation factor A protein-like 7 protein (1-100AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Transcription elongation factor A protein-like 7 (TCEA-like protein 7) (Transcription elongation factor S-II protein-like 7), TCEAL7
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	TCEAL7

## Image



### Western blot

All lanes: TCEAL7 antibody at 2µg/ml + A549 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 13 kDa

Observed band size: 13 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA863109LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA863109LA01HU at dilution of 1:100