




TCEAL7 Antibody

Product Code	CSB-PA863109LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRU2
Immunogen	Recombinant Human Transcription elongation factor A protein-like 7 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription elongation factor A protein-like 7 (TCEA-like protein 7) (Transcription elongation factor S-II protein-like 7), TCEAL7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TCEAL7

Image



Western blot

All lanes: TCEAL7 antibody at 2µg/ml + A549 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 13 kDa

Observed band size: 13 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA863109LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA863109LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.