



TCP10 Antibody, FITC conjugated

Product Code	CSB-PA614256HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12799
Immunogen	Recombinant Human T-complex protein 10A homolog protein (134-353AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-complex protein 10A homolog, TCP10, TCP10A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TCP10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.