





TGFB3 Antibody

Product Code	CSB-PA023449LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10600
Immunogen	Recombinant Human Transforming growth factor beta-3 proprotein protein (24-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-3 proprotein (TGF-beta-3) [Cleaved into: Latency-associated peptide (LAP)] (Transforming growth factor beta-3), TGFB3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TGFB3

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA023449LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA023449LA01HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.