



THTPA Antibody, Biotin conjugated

Product Code	CSB-PA883600LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BU02
Immunogen	Recombinant Human Thiamine-triphosphatase protein (2-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thiamine-triphosphatase (ThTPase) (EC 3.6.1.28), THTPA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	THTPA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.