



THY1 Antibody, Biotin conjugated

Product Code	CSB-PA023525LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04216
Immunogen	Recombinant Human Thy-1 membrane glycoprotein protein (20-130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thy-1 membrane glycoprotein (CDw90) (Thy-1 antigen) (CD antigen CD90), THY1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	THY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.