



TIFA Antibody, HRP conjugated

Product Code	CSB-PA839301LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CG3
Immunogen	Recombinant Human TRAF-interacting protein with FHA domain-containing protein A protein (1-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TRAF-interacting protein with FHA domain-containing protein A (Putative MAPK-activating protein PM14) (Putative NF-kappa-B-activating protein 20) (TRAF2-binding protein), TIFA, T2BP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TIFA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.