




TIMM8B Antibody

Product Code	CSB-PA897567LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5J9
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim8 B protein (33-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import inner membrane translocase subunit Tim8 B (DDP-like protein) (Deafness dystonia protein 2), TIMM8B, DDP2 DDPL TIM8B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TIMM8B

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA897567LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.