



TIMM8B Antibody, Biotin conjugated

Product Code	CSB-PA897567LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5J9
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim8 B protein (33-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import inner membrane translocase subunit Tim8 B (DDP-like protein) (Deafness dystonia protein 2), TIMM8B, DDP2 DDPL TIM8B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TIMM8B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.