





TIMM9 Antibody

Product Code	CSB-PA897296LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5J7
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim9 protein (2-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import inner membrane translocase subunit Tim9, TIMM9, TIM9 TIM9A TIMM9A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TIMM9

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA897296LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA897296LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.