



TIMM9 Antibody, Biotin conjugated

Product Code	CSB-PA897296LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5J7
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim9 protein (2-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import inner membrane translocase subunit Tim9, TIMM9, TIM9 TIM9A TIMM9A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TIMM9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.