



TMBIM1 Antibody, Biotin conjugated

Product Code	CSB-PA853376LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969X1
Immunogen	Recombinant Human Protein lifeguard 3 protein (1-109AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein lifeguard 3 (Protein RECS1 homolog) (Transmembrane BAX inhibitor motif-containing protein 1), TMBIM1, LFG3 RECS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMBIM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.