



TMED2 Antibody, HRP conjugated

Product Code	CSB-PA05049B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15363
Immunogen	Recombinant Human Transmembrane emp24 domain-containing protein 2 protein (26-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane emp24 domain-containing protein 2 (Membrane protein p24A) (p24) (p24 family protein beta-1) (p24beta1), TMED2, RNP24
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMED2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.