





TMEM115 Antibody

Product Code	CSB-PA023683LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12893
Immunogen	Recombinant Human Transmembrane protein 115 protein (187-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 115 (Placental protein 6) (Protein PL6), TMEM115, PL6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMEM115


Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA023683LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA023683LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA023683LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.