



TMEM14A Antibody, FITC conjugated

Product Code	CSB-PA023720LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6G1
Immunogen	Recombinant Human Transmembrane protein 14A protein (45-78AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 14A, TMEM14A, C6orf73
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM14A