



TMEM176B Antibody, FITC conjugated

Product Code	CSB-PA023758LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3YBM2
Immunogen	Recombinant Human Transmembrane protein 176B protein (1-64AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 176B (Protein LR8), TMEM176B, LR8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM176B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.