





TMEM217 Antibody

Product Code	CSB-PA023811LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N7C4
Immunogen	Recombinant Human Transmembrane protein 217 protein (150-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 217, TMEM217, C6orf128
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM217

Image



Western blot
 All lanes: TMEM217 antibody at 10µg/ml + Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 27, 25 kDa
 Observed band size: 27 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA023811LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.