



TMEM217 Antibody, HRP conjugated

Product Code	CSB-PA023811LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N7C4
Immunogen	Recombinant Human Transmembrane protein 217 protein (150-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 217, TMEM217, C6orf128
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM217