



# TMEM255B Antibody, HRP conjugated

Product Code	CSB-PA819893LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV15
Immunogen	Recombinant Human Transmembrane protein 255B protein (221-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 255B (Protein FAM70B), TMEM255B, FAM70B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM255B