



TMEM38B Antibody, HRP conjugated

Product Code	CSB-PA878900LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVV0
Immunogen	Recombinant Human Trimeric intracellular cation channel type B protein (249-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trimeric intracellular cation channel type B (TRIC-B) (TRICB) (Transmembrane protein 38B), TMEM38B, C9orf87
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TMEM38B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.