





TMEM45A Antibody

Product Code	CSB-PA023842LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWC5
Immunogen	Recombinant Human Transmembrane protein 45A protein (171-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 45A (DNA polymerase-transactivated protein 4) (Dermal papilla-derived protein 7), TMEM45A, DERP7 DNAPTP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM45A


Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA023842LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA023842LA01HU at dilution of 1:100



Immunofluorescent analysis of A2780 cells using CSB-PA023842LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.