



TMEM45A Antibody, Biotin conjugated

Product Code	CSB-PA023842LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWC5
Immunogen	Recombinant Human Transmembrane protein 45A protein (171-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 45A (DNA polymerase-transactivated protein 4) (Dermal papilla-derived protein 7), TMEM45A, DERP7 DNAPTP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM45A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.