




PIP4P2 Antibody

Product Code	CSB-PA023855LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4L2
Immunogen	Recombinant Human Type 2 phosphatidylinositol 4,5-bisphosphate 4-phosphatase protein (1-191AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PIP4P2

Image



Western Blot
 Positive WB detected in: Mouse Liver tissue lysate (30µg)
 All lanes: TMEM55A antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/20000 dilution
 Predicted band size: 29 kDa
 Observed band size: 29 kDa
 Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.