



# PIP4P2 Antibody, Biotin conjugated

Product Code	CSB-PA023855LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4L2
Immunogen	Recombinant Human Type 2 phosphatidylinositol 4,5-bisphosphate 4-phosphatase protein (1-191AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Type 2 phosphatidylinositol 4,5-bisphosphate 4-phosphatase (Type 2 PtdIns-4,5-P2 4-Ptase) (EC 3.1.3.78) (PtdIns-4,5-P2 4-Ptase II) (Transmembrane protein 55A), PIP4P2, TMEM55A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PIP4P2