



TMEM63A Antibody, Biotin conjugated

Product Code	CSB-PA023865LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94886
Immunogen	Recombinant Human CSC1-like protein 1 protein (212-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CSC1-like protein 1 (Transmembrane protein 63A), TMEM63A, KIAA0489 KIAA0792
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TMEM63A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.