



TMPRSS6 Antibody, HRP conjugated

Product Code	CSB-PA808531LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IU80
Immunogen	Recombinant Human Transmembrane protease serine 6 protein (98-301AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protease serine 6 (EC 3.4.21.-) (Matriptase-2), TMPRSS6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMPRSS6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.