



Tnfrsf4 Antibody, Biotin conjugated

Product Code	CSB-PA023981HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47741
Immunogen	Recombinant Mouse Tumor necrosis factor receptor superfamily member 4 protein (20-211AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 4 (OX40 antigen) (OX40L receptor) (CD antigen CD134), Tnfrsf4, Ox40 Txgp1
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Tnfrsf4