



TNFSF12 Antibody, Biotin conjugated

Product Code	CSB-PA023987LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43508
Immunogen	Recombinant Human Tumor necrosis factor ligand superfamily member 12 protein (43-149AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 12 (APO3 ligand) (TNF-related weak inducer of apoptosis) (TWEAK) [Cleaved into: Tumor necrosis factor ligand superfamily member 12, membrane form; Tumor necrosis factor ligand superfamily member 12, secreted form], TNFSF12, APO3L DR3LG
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFSF12