



TRAPPC3 Antibody, Biotin conjugated

Product Code	CSB-PA024229LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43617
Immunogen	Recombinant Human Trafficking protein particle complex subunit 3 protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trafficking protein particle complex subunit 3 (BET3 homolog), TRAPPC3, BET3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRAPPC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.