



Trem2 Antibody, Biotin conjugated

Product Code	CSB-PA024405LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99NH8
Immunogen	Recombinant Mouse Triggering receptor expressed on myeloid cells 2 protein (19-171AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Triggering receptor expressed on myeloid cells 2 (TREM-2) (Triggering receptor expressed on monocytes 2), Trem2, Trem2a Trem2b Trem2c
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Trem2