



Trem2 Antibody, Biotin conjugated

Product Code	CSB-PA024405LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99NH8
Immunogen	Recombinant Mouse Triggering receptor expressed on myeloid cells 2 protein (19-171AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Triggering receptor expressed on myeloid cells 2 (TREM-2) (Triggering receptor expressed on monocytes 2), Trem2, Trem2a Trem2b Trem2c
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Trem2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.