



TREML1 Antibody, HRP conjugated

Product Code	CSB-PA771477LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YW5
Immunogen	Recombinant Human Trem-like transcript 1 protein (16-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trem-like transcript 1 protein (TLT-1) (Triggering receptor expressed on myeloid cells-like protein 1), TREML1, TLT1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TREML1