



TRIM15 Antibody, Biotin conjugated

Product Code	CSB-PA887181LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C019
Immunogen	Recombinant Human Tripartite motif-containing protein 15 protein (122-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 15 (RING finger protein 93) (Zinc finger protein 178) (Zinc finger protein B7), TRIM15, RNF93 ZNF178 ZNFB7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.