



# TRIM15 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA887181LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9C019
<b>Immunogen</b>	Recombinant Human Tripartite motif-containing protein 15 protein (122-387AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Tripartite motif-containing protein 15 (RING finger protein 93) (Zinc finger protein 178) (Zinc finger protein B7), TRIM15, RNF93 ZNF178 ZNFB7
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	TRIM15