




TRIM17 Antibody

Product Code	CSB-PA897559LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y577
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM17 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRIM17

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, A549 whole cell lysate


All lanes: TRIM17 antibody at 1:500

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 55,40 kDa

Observed band size: 55 kDa



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA897559LA01HU at dilution of 1:67



Immunofluorescence staining of HepG2 cells with CSB-PA897559LA01HU at 1:25, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.