



TRIM24 Antibody, HRP conjugated

Product Code	CSB-PA024460LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15164
Immunogen	Recombinant Human Transcription intermediary factor 1-alpha protein (891-1012AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription intermediary factor 1-alpha (TIF1-alpha) (EC 2.3.2.27) (E3 ubiquitin-protein ligase TRIM24) (RING finger protein 82) (RING-type E3 ubiquitin transferase TIF1-alpha) (Tripartite motif-containing protein 24), TRIM24, RNF82 TIF1 TIF1A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.