



TRIM38 Antibody, FITC conjugated

Product Code	CSB-PA024475LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00635
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM38 protein (45-259AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM38 (EC 2.3.2.27) (RING finger protein 15) (RING-type E3 ubiquitin transferase TRIM38) (Tripartite motif-containing protein 38) (Zinc finger protein RoRet), TRIM38, RNF15 RORET
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRIM38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.