



TRIM51 Antibody, Biotin conjugated

Product Code	CSB-PA871567LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSJ1
Immunogen	Recombinant Human Tripartite motif-containing protein 51 protein (160-452AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 51 (SPRY domain-containing protein 5), TRIM51, SPRYD5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRIM51
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.