



TRIP13 Antibody, Biotin conjugated

Product Code	CSB-PA623008LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15645
Immunogen	Recombinant Human Pachytene checkpoint protein 2 homolog protein (181-390AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pachytene checkpoint protein 2 homolog (Human papillomavirus type 16 E1 protein-binding protein) (16E1-BP) (HPV16 E1 protein-binding protein) (Thyroid hormone receptor interactor 13) (Thyroid receptor-interacting protein 13) (TR-interacting protein 13) (TRIP-13), TRIP13, PCH2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TRIP13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.