




TRIT1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA867142LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H3H1 |
| Immunogen | Recombinant Human tRNA dimethylallyltransferase protein (1-324AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | tRNA dimethylallyltransferase (EC 2.5.1.75) (Isopentenyl-diphosphate:tRNA isopentenyltransferase) (IPP transferase) (IPPT) (hGRO1) (tRNA isopentenyltransferase) (IPTase), TRIT1, IPT MOD5 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | TRIT1 |

Image



Western blot

All lanes: TRIT1 antibody at 3µg/ml

Lane 1: Jurkat whole cell lysate


Lane 2: Raji whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 53, 38, 50, 44, 19 kDa

Observed band size: 53 kDa



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA867142LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.