



TRNAU1AP Antibody, HRP conjugated

Product Code	CSB-PA882127LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX07
Immunogen	Recombinant Human tRNA selenocysteine 1-associated protein 1 protein (33-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA selenocysteine 1-associated protein 1 (SECp43) (tRNA selenocysteine-associated protein 1), TRNAU1AP, SECP43 TRSPAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRNAU1AP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.