



# TRUB1 Antibody, Biotin conjugated

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA851978LD01HU                                                               |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q8WWH5                                                                           |
| Immunogen           | Recombinant Human Probable tRNA pseudouridine synthase 1 protein (2-349AA)       |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Human                                                                            |
| Tested Applications | ELISA                                                                            |
| Form                | Liquid                                                                           |
| Conjugate           | Biotin                                                                           |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified                                                         |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Alias               | Probable tRNA pseudouridine synthase 1 (EC 5.4.99.-), TRUB1, PUS4                |
| Immunogen Species   | Homo sapiens (Human)                                                             |
| Research Area       | Epigenetics and Nuclear Signaling                                                |
| Target Names        | TRUB1                                                                            |