



TRUB1 Antibody, HRP conjugated

Product Code	CSB-PA851978LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWH5
Immunogen	Recombinant Human Probable tRNA pseudouridine synthase 1 protein (2-349AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable tRNA pseudouridine synthase 1 (EC 5.4.99.-), TRUB1, PUS4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRUB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.