



TSACC Antibody, HRP conjugated

Product Code	CSB-PA842631LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A04
Immunogen	Recombinant Human TSSK6-activating co-chaperone protein (1-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TSSK6-activating co-chaperone protein (SSTK-interacting protein) (SSTK-IP), TSACC, C1orf182
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TSACC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.